

# CANINE UNDERWATER TREADMILL COMPARISON CHART

<b>There is No “NO” in Westcoast</b>		
	<b>Westcoast Treadmill</b>	<b>Competitive Treadmill</b>
Smooth Start-up, “No jerk” for the dogs	Yes	No
Speeds from a fraction to 8 mph	Yes	No- .5 to 5 mph
Less Expensive electrical pull	Yes	No
Mobility – Entire system in on locking wheels.	Yes	No - an additional installation charge to move the equipment
Requires zero water pressure to ‘hold’ a water tight seal	Yes	No - requires a level of water to hold the door shut
Runs Quietly	Yes	No
Fits though normal doorways (36” wide)	Yes	No - requires on-site installation within the treatment room
Requires No special wiring. 220 VAC, and 30 amp circuit. Same as clothes dryer	Yes	No - requires 3 220 VAC dedicated circuits Avg. cost \$3000
Door Swings Out – Can accommodate “longer” dogs	Yes	No
Stainless steel, easy to clean	Yes	No
Sloped sides next to belt	Yes	No - Flat - Dogs can cheat often mandating that a technician be in the unit at all times
Pumps can be up to 75 feet away	Yes	No
Glass sided for viewing	Yes	No
Glass sided for easy cleaning	Yes	No
Glass sided for scratch resistance	Yes	No
No special tools for installation	Yes	No
Unit plugs into one outlet	Yes	No
Sling Support Bars for debilitated animals	Yes	No
Sand filtered water	Yes	No
Water Height Gauge on Side	Yes	No
Less then \$1 per day of electricity use	Yes	No
Can easily clean under the belt	Yes	No – Requires wrenches and excessive down time
Can dispose of water in use, after defecation.	Yes	No
Can immediately use treadmill after defecation.	Yes	No - Entire system must be refilled & heated
Quiet Pumps, Filters and Belt	Yes	No



**Joseph P. Thompson & Associates, LLC**

**866-646-8684**