

# PiezoWave<sup>2</sup> VET

## SMALL ANIMAL SHOCKWAVE THERAPY COMPARISON CHART

(for F7 G3 Therapy Source)	PiezoWave Vet	VersaTron 4 Paws
Multiple Therapy Heads (Trodes) Required	NO	YES
Therapy Depth of Therapy Ranges	8 - 0mm to 30mm in 5mm increments	2 - 5mm and 20mm
Linear Therapy Source For Wounds Available	YES	NO
Complete Animal Sedation Required	NO	YES
Ear Protection Required	NO	YES
Unit Cost	\$30,000	\$30,000
Cost Per Shockwave	.0013 cents	9.5 cents
Cost Per 1000 Shocks	\$1.30	\$95.00
Shocks Per Therapy Source	5 Million +	20,000
Number of Therapy Source Units/5 Million Shocks	1	250
Cost of Replacement Therapy Source	\$6,500	\$1,900
<b>Total Cost of Ownership Based on 5M Shocks</b>	<b>\$30,000</b>	<b>\$491,200</b>



[www.mtavet.com](http://www.mtavet.com)

866.646.8684